T/PUMP® SERIES LOCALIZED TEMPERATURE THERAPY SYSTEM

REF TP700

REF TP700C

Symbols


	Operating instructions/Consult instructions for use
	General warning
	Caution
	Do not puncture
	Warning non-ionizing radiation
	Catalogue number
	Serial number
	Manufacturer
	For US Patents see www.C2Dx.com/patents
	Dangerous voltage
	Type BF applied part
	Medical Equipment Recognized by Underwriters Laboratories LLC With Respect to Electric Shock, Fire, and Mechanical Hazards only in accordance with UL 60601-1: 2003 and CAN/CSA-C22.2 No. 601-1, ASTM F 2196-2002
	Warm water fill line
	Cold water fill line
	Water levels
	Water flow
	 CAUTION Federal law (USA) restricts this device to sale by or on the order of a physician.

Table of Contents

Warning / Caution / Note Definition	2
Summary of safety precautions.....	3
Introduction	4
Product description.....	4
Indications for use.....	4
Expected service life.....	4
Contraindications.....	5
Specifications.....	6
Product illustration.....	7
Product keypad.....	8
Contact information.....	8
Serial number location.....	9
Operation.....	10
Connecting single or multiple pads.....	10
Connecting and disconnecting the Clik-Tite® Connectors.....	10
Connecting and disconnecting the Colder Connectors.....	11
Opening and closing hose pinch clamps.....	11
Starting the pump.....	12
Starting therapy cycles.....	13
Stopping the pump.....	13
Storing the pump, short term.....	14
Storing the pump, long term.....	14
Draining the pump.....	14
Cleaning.....	15
Cleaning the external surface.....	15
Cleaning the water circuit and hoses.....	15
Preventive maintenance	16
Accessories.....	17
EMC Information.....	19
Warranty	23
Warranty exclusion and damage limitations.....	23
Return authorization.....	23
Damaged product.....	23
International warranty clause.....	23

Warning / Caution / Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note: Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Carefully read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Always fill the pump reservoir with water before operating. Do not plug in the product until you fill the water reservoir.
 - Always plug the product into a hospital grade wall outlet.
 - Do not attach the power cord to any moving parts of the bed frame or any location that may pinch or cause power cord damage. Improper routing can cause power cord damage and possible electrical shock.
 - Always set the temperature as directed by a physician.
 - Always use the **T/Pump® Temperature Therapy Pump** with **Mul-T-Pad® Temperature Therapy Pad** accessories. Do not use with competitor pads.
 - Risk of explosion. Do not use in the presence of flammable anesthetics.
 - Always check the skin integrity of the body surface that therapy is applied. Evaluate patient response to the temperature application.
 - Always check the patient's skin for adverse reactions every 30 minutes or as directed by the physician.
 - Do not place additional heat sources between the patient and the pad. Skin damage may result.
 - Always allow for additional monitoring of the following groups or conditions. Pediatric patients. Patients with impaired circulation (Ischemia). Closely monitor the areas of body that are beneath the thermal transfer devices (Ischemia).
 - Always unplug the product before cleaning or disinfecting to avoid the risk of shock.
 - Risk of electric shock. Disconnect power before servicing or cleaning the **T/Pump® Temperature Therapy Pump**.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Federal law (USA) restricts this device to sale by or on the order of a physician.
 - Do not cover the product with blankets, pillows, or other insulating materials. Air flow is required to maintain system performance.
 - Closely monitor the areas of body that are beneath the thermal transfer devices. Also check for areas of the body for use of topical solutions, or other heat sources. These may interfere with the treatment prescribed.
 - Always set the fluid temperature and monitor the temperature of the patient and patient skin condition as directed by a physician.
 - Do not use quaternary disinfectant solutions on the external surface of this product.
 - Do not spray cleaning or disinfectant solutions onto the pump. Frequent or prolonged disinfectant solutions may result in damage to the product.
 - Do not use the product next to or stacked with other equipment. If next to or stacked is necessary, observe the product to verify normal operation in the configuration in which it used.
-

Introduction

This manual assists you with the operation or maintenance of your C2Dx® product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note:

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- C2Dx® continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact C2Dx® Customer Service or Technical Support at 888-902-2239.

Product description

The **T/Pump®** Localized Temperature Therapy system delivers temperature controlled water. The water flows through an impeller driven pump through a bi-lumen hose. The hose connects to thermal transfer devices that can be applied over, under, or around the patient. Water temperature increases with the use of the heater. Water temperature decreases with the use of ice. The system delivers and removes thermal energy to and from the local anatomical sites. This occurs at the contact point of the patient skin and thermal transfer device to raise or lower the temperature of local anatomical sites.

Indications for use

The **T/Pump®** localized temperature therapy system is for use in supplying warm or cold water at controlled temperatures via thermal transfer devices. This is for the application of temperature therapy in situations where a physician determines that temperature therapy is necessary. Intended for adult and pediatric human patient's skin as attended by health care professionals and lay persons. For use in healthcare settings and home use.

Localized temperature therapy is of particular benefit in treating the following applications:

- Orthopedic conditions such as acute injuries, chronic pain, and muscle spasm and strains.
- Skin trauma such as abscesses, boils, bruises, burns, and contusions.
- Other medical conditions such as chronic arthritis, neuritis, phlebitis, tendonitis and IV infiltration, and symptoms such as infections and localized pain.

For use with non-acute traumatized tissue, impaired mental status, insensate body surface, and oxygen therapy. If an oxygen tent is in use, do not use the pump inside the tent.

Note: Applications are as determined by a physician.

Expected service life

The **T/Pump® Temperature Therapy Pump** controller has a two year expected service life under normal use conditions and with appropriate periodic maintenance.

Introduction

Contraindications

C2Dx® promotes the assessment of each patient and appropriate usage by the operator. Contraindications for this product are for patients with the conditions specified below.

Heating contraindications include:

- Application to a body surface with compromised blood flow (Ischemia, area under pressure, arterial insufficiency)
- Application to a patient with an increased tendency to bleed (aggravates potential for hemorrhage)
- Application to a body surface with possible malignancy (tissue metabolism is increased and therefore, the growth potential of the malignant tissues)
- Treatment of hematoma within the first 24 to 48 hours (potential for re-bleeding and hemorrhage). Recent sprain or fracture (acute inflammatory response).
- In combination with topical solutions whose toxicity may be affected by the application of heat
- In combination with other heat sources



Cooling contraindications include:

- Application to a body surface with compromised blood flow (Ischemia, area under pressure, arterial insufficiency)
- Application to body surface with known vascular impairment such as frostbite, arteriosclerosis or ischemia
- Application to body surface in patients with hypersensitivity to cold, such as people with Raynaud's phenomenon, cold urticaria, cryoglobulinemia, and paroxysmal cold hemoglobinuria
- Application to body surface in patients with impaired sensation
- In combination with topical solutions whose toxicity may be affected by the application of cold

Introduction

Specifications


Models	TP700 / TP700C	
Temperature setpoints Fahrenheit	107° F, 100° F, 95° F, 50° F	
Temperature setpoints Celsius	42° C, 38° C, 35° C, 10° C	
Therapy duration	20 or 30 minutes or continuous cycle	
Classification	Class I equipment with Type BF applied part suitable for continuous operation. Not classified for protection against ingress of liquid. Not classified for use in the presence of flammable anesthetics.	
Size	11.5 in. x 8 in. x 8 in.	29.2 cm x 20.3 cm x 20.3 cm
Weight	Reservoir empty	6.5 lb 2.9 kg
	Reservoir filled with water to the heating level	9 lb 4.0 kg
Maximum reservoir capacity	93 oz	2.75 l
Flow rate, minimum with pad attached	9 gph	34 lph
Average temperature accuracy	±2° F at 107° F	±1° C to 42° C
Maximum contact surface temperature	107° F	42° C
High limit safety temperature	110° F to 120° F	43.3° C to 49° C
Power cord	International (harmonized) 3 wire cord set	
Dual hose	10 ft	305 cm
Current leakage, maximum	100 µA	
Ground resistance, maximum	0.5 mΩ	
Electrical requirements	Voltage	120 VAC ±10%
	Frequency:	60 Hz
	Current in amperes	3.1

Environmental conditions	Operation	Storage and transportation
Ambient operating temperature	 <p>60°F (15.6°C) to 90°F (32.2°C)</p>	 <p>-20°F (-28°C) to 120°F (48°C)</p>

C2Dx® reserves the right to change specifications without notice.

Introduction

Product illustration


A	Handle vents
B	Keypad
C	Flow indicator
D	Warm or cool delivery
E	Therapy cycle times
F	Hose connections

G	Hose storage
H	Cord storage
I	Tethered open cap
J	Handle
K	On/Standby button

Introduction

Product keypad

A	Indicator / Warning light	Red warning LED flashes with a short audible beep to inform the operator that the product has reached the desired temperature or is moving to the next cycle
B	Water flow	Check hoses or clamps for kinks or occlusions
C	Water level	Check water level
D	Setpoints	Toggle through the four setpoints, temperature in C and F
E	On/Standby button	Green indicates the unit is on / Yellow indicates power is supplied to the unit but the unit is not on
F	Therapy cycles	Continuous cycle, 30 minute cycle, or 20 minute cycle
G	Setpoint lock	Prevents tampering of settings

Contact information

Contact C2Dx® Customer Service or Technical Support at 888-902-2239.

C2Dx® Inc
 555 E Eliza St, Ste A
 Schoolcraft, MI 49087
 USA

To view your operations or maintenance manual online, see www.C2Dx.com/manuals

Have the serial number (A) of your C2Dx® product available when calling C2Dx® Customer Service or Technical Support. Include the serial number in all written communication.

Introduction

Serial number location

Operation

Connecting single or multiple pads

The **Mul-T-Pad®** products provide the interface for delivering temperature therapy. The unique button design allows water to flow and provides trouble-free operation when the pad is folded to form a customized fit. You apply the pads to the part of the body that require therapy. The circulating water maintains the pad at the setpoint temperature.

The pad is connected to the pump with **Clik-Tite® Connectors** or Colder connectors.

You can interconnect the **Mul-T-Pad® Temperature Therapy Pad** (Figure 1-1 on page 10) by using the **Clik-Tite® Connectors** to provide therapy to more than one body part at a time (Figure 1-2 on page 10).

Note: Colder style connectors require Adaptor P/N 77926004 to interconnect multiple pads.


Figure 1-1: Localized temperature system with single pad


Figure 1-2: Localized temperature therapy system with multiple pads

Connecting and disconnecting the Clik-Tite® Connectors

TP700 has **Clik-Tite®** style connectors.

To connect the **Clik-Tite® Connectors** from the hose to the pad:

1. Insert the male fittings into the female fittings with a twisting motion (Figure 1-3 on page 10 and Figure 1-4 on page 10).


Figure 1-3: Connecting fittings


Figure 1-4: Twist fittings

2. After you insert fittings (Figure 1-5 on page 10), snap the lock ring into place (Figure 1-6 on page 10)


Figure 1-5: Insert fittings


Figure 1-6: Lock ring

Operation

Connecting and disconnecting the Clik-Tite® Connectors (Continued)

3. To disconnect, reverse the procedure.

Connecting and disconnecting the Colder Connectors

TP700C has **Colder** style connectors

To connect Colder connectors from the hose to the pad:

1. Insert the male coupling into the female coupling (Figure 1-7 on page 11). When you hear an audible click, you have joined the connectors.

Figure 1-7: Connecting couplings

2. To disconnect, press down on the thumb tab (A) of the female coupling. The couplings will partially disconnect.
3. Pull the male coupling out to disconnect (Figure 1-8 on page 11).


Figure 1-8: Disconnecting couplings

Opening and closing hose pinch clamps

To open or close the hose pinch clamps:

1. Open the clamp by pushing on the serrated end (Figure 1-9 on page 11).


Figure 1-9: Open the clamp


Figure 1-10: Close the clamp

2. Close the clamp by pressing the clamp together (Figure 1-10 on page 11).

Operation

Starting the pump

WARNING

- Always fill the pump reservoir with water before operating. Do not plug in the product until you fill the water reservoir.
 - Always plug the product into a hospital grade wall outlet.
 - Do not attach the power cord to any moving parts of the bed frame or any location that may pinch or cause power cord damage. Improper routing can cause power cord damage and possible electrical shock.
 - Always set the temperature as directed by a physician.
 - Always use the **T/Pump® Temperature Therapy Pump** with **Mul-T-Pad® Temperature Therapy Pad** accessories. Do not use with competitor pads.
-

CAUTION

- Federal law (USA) restricts this device to sale by or on the order of a physician.
 - Do not cover the product with blankets, pillows, or other insulating materials. Air flow is required to maintain system performance.
 - Closely monitor the areas of body that are beneath the thermal transfer devices. Also check for areas of the body for use of topical solutions, or other heat sources. These may interfere with the treatment prescribed.
-

To start the pump:

1. Before filling the pump reservoir, attach a pad to the connector hose (see [Connecting single or multiple pads on page 10](#)). Make sure that there are no kinks in the hose or pad.
2. Open the hose clamps.
3. Open the fill cap on the top of the pump.
4. To fill for cooling
 - Fill with cold water to the Cooling water line.
 - Fill with ice to the full capacity of the reservoir.

Note: When using the pump for cooling, the ice will eventually melt and the setpoint will start flashing.

5. To fill for heating, fill with room temperature water to the heating water line.
6. Plug the pump into a properly grounded hospital grade wall outlet.
7. Press the On/Standby button. The light next to the selected temperature begins to flash.
8. Use the keypad to set the temperature as directed by the physician.

Notes

- After setting the temperature, press and hold the Setpoint lock button for two seconds to lock the setpoint.
 - If you go past the desired setpoint, keep pressing the button to start at the beginning of the setpoint values.
 - It takes about 15 minutes to reach the selected water temperature. The light next to the selected temperature becomes steady.
 - Check the water level. If the water level drops below the operating level, add water.
9. Apply the **Mul-T-Pad® Temperature Therapy Pad** to the patient as directed by a physician.
 10. Position the pump at or above the level of the pad.

Note: Do not place the pump below the pads, water will drain into the pump during shut off. If the pump has been overfilled or if you connected multiple pads, excess water can leak.

Operation

Starting therapy cycles

WARNING

- Risk of explosion. Do not use in the presence of flammable anesthetics.
 - Always check the skin integrity of the body surface that therapy is applied. Evaluate patient response to the temperature application.
 - Always check the patient's skin for adverse reactions every 30 minutes or as directed by the physician.
 - Do not place additional heat sources between the patient and the pad. Skin damage may result.
 - Always allow for additional monitoring of the following groups or conditions. Pediatric patients. Patients with impaired circulation (Ischemia). Closely monitor the areas of body that are beneath the thermal transfer devices (Ischemia).
-

CAUTION

- Always set the fluid temperature and monitor the temperature of the patient and patient skin condition as directed by a physician.
-

Heat therapy is effective in the dilation of blood vessels, which increases the blood flow to the heated area.

Cooling therapy assists in the vasoconstriction, decreasing blood flow, and decreasing the metabolism in the affected area.

To start therapy cycles

1. Press the On/Standby button to turn the pump on.
2. Set the temperature setpoint. The unit will start to heat up.
3. The off time period is equal to the Therapy Cycle selected. Continuous does not have an off time period.
4. Select the warming period for a 20 minute or 30 minute cycle. The LED turns solid and the selected temperature LED flashes.:

Notes

- The red warming LED flashes with a short audible beep to inform the operator the unit is at the set temperature. The timed therapy is starting.
- When the 20 minute or 30 minute time period is completed, the heater and pump shuts off. The water flow to the pads is stopped. The setpoint and therapy time LEDs blink to indicate the off time period.
- The off time period is equal to the Therapy Cycle selected. Continuous does not have an off time period
- After the off time period is completed, the pump and heater restart.
- The time therapy LED goes solid and the setpoint LED flashes.
- After the setpoint is reached, the red LED flashes with a short audible beep to let the operator know the next therapy cycle has started.

Stopping the pump

To stop the pump:

1. Press the On/Standby button. The standby light illuminates.
2. Unplug the product.
3. Close all of the hose clamps.
4. Disconnect the pads from the pump.

Note:

- There will be a series of audible beeps if device loses power before On/Standby button is pressed.
- To prevent water spillage, raise the connections above the level of the pad and the pump before disconnecting the pad from the pump.

Operation

5. Coil the hose and attach the Clik-Tite® Connectors together on the hose.
6. Secure the hose to the pump by using the tube set strap.
7. Wrap the power cord around the product.

Storing the pump, short term

Short term storage is less than one day.

To store the pump short term:

1. Close the hose clamps.
2. Disconnect the pads.
3. Connect the ends of the connector hoses together, where applicable.
4. Open the hose clamps.
5. Leave water in the reservoir.
6. Coil and fasten the hose using the tube set strap and wrap the power cord around the product.

Note: Change the water monthly or more often depending on use.

Storing the pump, long term

Long term storage is defined as more than one day.

To store the pump long term:

1. Drain the pump ([see Draining the pump on page 14](#)).
2. Coil and fasten the hose using the tube set strap and wrap the power cord around the product.

Draining the pump

To drain the pump:

1. Disconnect the pump from AC power.
2. Clamp the hose clamps.
3. Disconnect the pads or hoses from one another, keep the hoses at or above the level of the pump.
4. Open the hose clamps.
5. Remove the fill cap and invert the pump over a sink to drain.
6. Replace the fill cap after all of the fluid has drained from the hoses and reservoir.
7. Connect the ends of the connector hoses together, where applicable.

Cleaning

Empty and clean the pump on a monthly basis following these procedures. Complete this process more often depending on use.

WARNING

Always unplug the product before cleaning or disinfecting to avoid the risk of shock.

CAUTION

- Do not use quaternary disinfectant solutions on the external surface of this product.
 - Do not spray cleaning or disinfectant solutions onto the pump. Frequent or prolonged disinfectant solutions may result in damage to the product.
-

Cleaning the external surface

Recommended cleaning agents include:

- A clean damp cloth and mild soapy water
- Phenolic disinfectant or 10% bleach solution

To clean the external side of the pump:

1. Unplug the controller.
2. Wipe the outside of the product with the following:
 - a. Using a clean, soft cloth moistened with a mild soap and water solution to remove foreign material
 - b. Using a clean, soft cloth moistened with a phenolic disinfectant or 10% bleach solution.
3. Follow specified contact time in accordance with the manufacturer's instructions for use.
4. Using a clean, soft cloth, wipe the outside of the product with water.
5. Dry the external surface with a clean, dry cloth.

Cleaning the water circuit and hoses

Tools Required:

- Disinfectant Cleaner 2.0 (8001-000-333)
 - Personal protective equipment, based on cleaning agent manufacturer's instructions
 - Soft, lint free cloth (2 or more)
1. Prepare a germicidal solution according to the manufacturer's instructions for the Disinfectant Cleaner 2.0 .
 2. Drain the pump.
 3. Connect the hose set together, where applicable, or attach a pad to the connector hose.
 4. Fill the reservoir to the heating water line on the back of the reservoir.
 5. Select 95° F (35° C) temperature on the keypad.
 6. Press the **On/Standby** button to start the pump.
 7. Circulate the solution for one hour.
 8. Drain the solution from the reservoir.
 9. Rinse and drain the reservoir with water.
 10. Dry the reservoir inside and outside with a dry, lint free cloth.

Preventive maintenance

Remove product from service before you perform the preventive maintenance inspection. Check all items listed during annual preventive maintenance for all C2Dx® Medical products. You may need to perform preventive maintenance checks more often based on your level of product usage. Service only by qualified personnel

WARNING

Always unplug the product before cleaning or disinfecting to avoid the risk of shock.

Note: Change the water monthly or more often depending on use.

Inspect the following items:

- _____ Line cord, the plug, and the housing for cuts or cracked insulation and quality of the strain relief, if damaged replace
- _____ Membrane panel for cuts or cracks, if damaged replace
- _____ Plastic parts on the pump for cracks
- _____ Connectors for cracks and missing O-rings, replace connectors if necessary
- _____ Hose connections and connectors function properly
- _____ Pads and hoses are free from tears, cracks, and water leaks

Refer to the **T/Pump® Temperature Therapy Pump** Maintenance Manual Functional Checks section that and complete during preventive maintenance.

Inspection	Value	Okay	Action Needed	Action Taken
Measure current leakage <100 µA				

Product serial number:
Completed by:
Date:

Accessories

These accessories may be available for use with your product. Confirm availability for your configuration or region. Call C2Dx® Customer Service: 888-902-2239.

Test tools

Catalog Number	Product Name	Number in Carton
TPT9	Measures the temperature and water flow	1
100925000	Test probe	1

Cleaner

Catalog Number	Product Name	Number in Carton
8001-000-333	Disinfectant cleaner 2.0	1

Optional accessories

Catalog Number	Product Name	Number in Carton
TPS2	Temperature therapy stand with storage basket and cord hooks	1
TPIV	IV pole mount plate	1

Notes

- Only use **Mul-T-Pad®** products with the **T/Pump® Temperature Therapy Pump**. The unique button design allows optimal water flow and provides trouble-free operation when the pad is folded.
- The **Mul-T-Pad®** series products with **Clik-Tite® Connectors** can be interconnected to provide therapy to more than one body site at a time.

Mul-T-Pad® Model	Connector	Length	Width	Expected life in days / single patient use	Material	Quantity per carton
8002-062-012	Clik-Tite® Connectors For use with TP700	18 in. (46 cm)	13 in. (33 cm)	30	Nonwoven fabric on one side, pliable polymer on the other side	20
8002-062-022	Clik-Tite® Connectors For use with TP700	22 in. (56 cm)	15 in. (38 cm)	30	Nonwoven fabric on one side, pliable polymer on the other side	20
8002-062-026	Clik-Tite® Connectors For use with TP700	26 in. (66 cm)	18 in. (46 cm)	30	Nonwoven fabric on one side, pliable polymer on the other side	10

Accessories

8002-062-222	Clik-Tite® Connectors For use with TP700	22 in. (56 cm)	15 in. (38 cm)	90 days / reusable	Heavy polymer, reusable	10
8002-062-612	Colder connectors For use with TP700C	18 in. (46 cm)	13 in. (33 cm)	30	Non-woven fabric on both sides	10
8002-062-622	Colder connectors For use with TP700C	22 in. (56 cm)	15 in. (38 cm)	30	Non-woven fabric on both sides	10
8002-062-626	Colder connectors For use with TP700C	26 in. (66 cm)	18 in. (46 cm)	30	Non-woven fabric on both sides	10

EMC Information

The **T/Pump® Temperature Therapy Pump** with EN60601-2: 2001 Second Edition

(CISPR Classified as Class B, Group 1 ISM equipment)

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following pages.

Portable and mobile RF communications equipment can effect Medical Electrical Equipment.

CAUTION

Do not use the product next to or stacked with other equipment. If next to or stacked is necessary, observe the product to verify normal operation in the configuration in which it used.

Guidance and Manufacturer's declaration - Electromagnetic Immunity

The **T/Pump®** system is suitable for use in the electromagnetic environment specified below. The customer or the user of the **T/Pump®** system should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality is that of typical commercial or hospital environment.
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11	$< 5\% U_t$ ($> 95\%$ dip in U_t) for 0,5 cycle $40\% U_t$ (60% dip in U_t) for 5 cycles $70\% U_t$ (30% dip in U_t) for 25 cycles $< 5\% U_t$ ($> 95\%$ dip in U_t) for 5 sec	$< 5\% U_t$ ($> 95\%$ dip in U_t) for 0,5 cycle $40\% U_t$ (60% dip in U_t) for 5 cycles $70\% U_t$ (30% dip in U_t) for 25 cycles $< 5\% U_t$ ($> 95\%$ dip in U_t) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the T/Pump® system requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.

EMC Information

(Continued)

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
--	-------	-------	---

Note:

UT is the a.c. mains voltage prior to applications of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the T/Pump®.

The **T/Pump®** system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **T/Pump®** system can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **T/Pump®** system as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1,2\sqrt{P}$	80 MHz to 800 MHz $d=1,2\sqrt{P}$	800 MHz to 2,5 GHz $d=2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note:

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note:

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

EMC Information

The **T/Pump**® system is suited for use in the electromagnetic environment specified below. The customer or the user of the **T/Pump**® system should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
<p>Conducted RF IEC 61000-4-6 *</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>3 Vrms 3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the T/Pump® system, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p style="text-align: center;">Recommended Separation Distance</p> <p style="text-align: center;"> $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}$ 80 MHz to 800 MHz $d=2,3\sqrt{P}$ 800 MHz to 2,5 GHz </p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

Note:
At 80 MHz and 800 MHz, the higher frequency range applies.

Note:
These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **T/Pump**® system is used exceeds the applicable RF compliance level above, the **T/Pump**® system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **T/Pump**® system.

^bOver the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.

EMC Information

Guidance and Manufacturer's declaration - Electromagnetic Emissions

The **T/Pump**® system is intended for use in an electromagnetic environment specified below. The customer or the user of the **T/Pump**® system should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment
RF Emissions CISPR 11	Group 1	The T/Pump ® system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The T/Pump ® system is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2 *	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

Warranty

C2Dx® Medical Division, a division of C2Dx® Corporation, warrants to the original purchaser the C2Dx® Model TP700 and TP700C **T/Pump® Temperature Therapy Pump**, to be free from defects in material and workmanship for a period of one year after date of delivery.

The above noted warranty period applies only to the original purchaser of the **T/Pump® Temperature Therapy Pump** product and begins on the date of delivery to such original purchaser.

Warranty exclusion and damage limitations

The express warranty set forth herein is the only warranty applicable to the product. Any and all other warranties, whether express or implied, including any implied warranty of merchantability or fitness for a particular purpose are expressly excluded by C2Dx®. In no event shall C2Dx® be liable for incidental or consequential damages.

Return authorization

Product cannot be returned without prior approval from the C2Dx® Customer Service Department. An authorization number will be provided which must be printed on the returned product. C2Dx® reserves the right to charge shipping and restocking fees on returned product. Special, modified, or discontinued products are not subject to return.

Damaged product

ICC Regulations require that claims for damaged product must be made within fifteen (15) days of receipt of the product. Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt. Upon prompt notification, C2Dx® will file a freight claim with the appropriate carrier for damages incurred. Claims will be limited in amount to the actual replacement cost. In the event that this information is not received by C2Dx® within the fifteen (15) day period following the delivery of the product, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full within thirty (30) days of receipt. Claims for any incomplete shipments must be made within thirty (30) days of invoice.

International warranty clause

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Contact your local C2Dx® Medical representative for additional information.

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