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
Figure 1: Batch lot code
Figure 1: Código de lote
Figure 1 : Code de lot


Figure 2: Connect one pad
Figure 2 : Conecte una almohadilla
Figure 2 : Connecter un coussinet

Figure 3: Connect multiple pads
Figure 3: Conecte varias almohadillas
Figure 3 : Connecter plusieurs coussinets

Introduction

This manual assists you with the operation of the C2Dx Model 8002-062-012, 8002-062-022, and 8002-062-026 Mul-T-Pad®. Read this manual before operating this product and keep a copy on file. Set methods and procedures to educate and train your staff on the safe operation of this product.

The T/Pump® localized therapy system supplies warm or cold water at controlled temperatures. Delivery of the water is through the thermal transfer devices (Mul-T-Pad®) for localized temperature therapy. The T/Pump is for use in situations where physicians determine heat therapy is necessary or desired.

Localized temperature therapy is of particular benefit in treating the following:

- Orthopedic conditions such as acute injuries, chronic pain, lower back pain, muscle spasms, and strains;
- Skin trauma such as abscesses, boils, bruises, burns and contusions;
- Other medical conditions such as chronic arthritis, neuritis, phlebitis, tendonitis, and IV infiltration; and symptoms such as infections and localized pain.

Product description

The Mul-T-Pad is a lightweight polymer with non-woven fabric on one side. The Click-Tite® connectors provide the means to connect to the temperature controller. The button design allows water to flow when the pad is folded to form a customized fit.

Expected life

The Mul-T-Pad has a 30 day single patient use expected life upon first use and under normal use conditions.

Contraindications

- For heating:
- Application to a body surface with compromised blood flow (ischemia, area under pressure, arterial insufficiency).
 - Application to a patient with an increased tendency to bleeding (aggravates potential for hemorrhage).
 - Application to a body surface with possibility of malignancy (tissue metabolism is increased and therefore, the growth potential of the malignant tissues).
 - Treatment of hematoma within first 24-48 hours (potential for re-bleeding and hemorrhage). Recent sprain or fracture (acute inflammatory response).
 - In combination with topical solutions whose toxicity may be affected by the application of heat.
 - In combination with other heat sources.

For cooling:

- Application to body surface with compromised blood flow (ischemia, area under pressure, arterial insufficiency)
- Application to body surface with known vascular impairment such as frostbite, arteriosclerosis or ischemia
- Application to body surface in people with hypersensitivity to cold, such as people with Raynaud's phenomenon, cold urticaria, cryoglobulinemia, and paroxysmal cold hemoglobinuria
- Application to body surface in people with impaired sensation
- In combination with topical solutions whose toxicity may be affected by the application of cold

Specifications

| | | | | | |
|------------------------|--|--------------|--------------|--------------|--------------|
| Material | rayon, ethylene vinyl acetate (EVA), polypropylene, polyethylene | | | | |
| Model | 8002-062-012 | 8002-062-022 | 8002-062-022 | 8002-062-026 | |
| Length | 18 in. | 45 cm | 22 in. | 56 cm | 26 in. 66 cm |
| Width | 13 in. | 33 cm | 15 in. | 38 cm | 18 in. 45 cm |
| Compatible controllers | T/Pump series temperature controllers | | | | |

| | | |
|--------------------------|--------------------------------|-------------------------------|
| Environmental Conditions | Operating | Storage and transportation |
| Temperature | 60°F (15.6°C) to 90°F (32.2°C) | -20°F (-29°C) to 120°F (48°C) |
| Relative humidity | 30% to 75% | 10% to 95% |

C2Dx reserves the right to change specifications without notice.

Symbols

- General warning
- Caution
- Consult instructions for use
- Do not puncture
- Quantity
- Single patient use
- Batch code
- Catalogue number
- Date of manufacture
- Manufacturer

Contact information

Contact C2Dx Customer Service or Technical Support at 1-888-902-2239.

C2Dx, Inc.
555 E. Eliza St. | Ste. A
Schoolcraft, MI 49087
USA

Have the Lot batch code (A) (Figure 1) available when you call Customer Service. Include the Lot batch code in all written communication.

Safety Precautions

WARNING

- Do not use if the package is damaged.
- Always use Mul-T-Pad products with an approved temperature controller. Using this product with other controllers has not been tested.
- Do not reuse this product on another patient to avoid the risk of cross-contamination and infection.
- Do not use this product to position or transfer a patient.
- Do not allow hoses to fold, kink, or to wrap around a patient.
- Always pre-fill the pad before applying to the patient.
- Always check the patient's skin condition of areas in contact with the pad at a minimum of 30 minutes or as directed by a physician.
- Always consult the controller operations manual before connecting this product.

CAUTION

- US Federal Law restricts this device to sale or by on the order of a licensed physician.
- Do not use sharp objects or pins with this product.
- Do not clean this product. If soiled, dispose of product according to hospital protocol. This is a single patient use product.
- Always inspect the product for tears, cuts, holes, stains, or any other damage before use. Always monitor more frequently pediatric patients and patients with impaired circulation.
- Always use minimal layers of sheeting and incontinence pads. Too many layers between the patient's skin and the pad will reduce the cooling or warming capabilities of the system.

Operations

Assessing the patient's skin condition

Note: Follow your facilities skin care protocol. The following instructions are recommendations.

1. Document initial skin assessment before therapy.
2. Check for any foreign items such as medical patches, IV's, or sutures.
3. Check for fluids such as moisture or gels.
4. Verify the pad size is correct and pre-filled before applying to the patient.

Note: Apply the pad directly to the skin or add minimal layers.

Connecting and pre-filling the pad (Click-Tite connectors)

1. Close the pinch clamps on the connector hose and pad (Figure 4).
2. Insert the male coupling of pad into the female coupling of hose with a twisting motion. Press until you hear a click (Figure 5 and 6).
3. Snap (Figure 7) and Lock (Figure 8) the fitting.
4. With a light pull, make sure that you locked the connectors.
5. Open all clamps on the connector hose and the pad (Figure 9).
6. Insert the power cord into a wall outlet.
7. Press the power button to turn on the controller. See the controller instructions for use.
8. Check that the pad fills with water.

Notes:

- Iodine based disinfectants will stain the product. Staining does not affect the patient nor the use of the product.
- Do not use pins or sharp objects with this product.
- If you fold the pad, place the end of the pad that contains the tubing away from the patient.
- Interconnect using the Click-Tite® connectors to provide therapy to more than one body site at a time (Figure 2 and Figure 3).

Rechecking the patient's skin condition

Recheck the entire area of the patient's skin condition that is in contact with the pad at least every 30 minutes. Note any change in the skin integrity that relates to:

- Excessive moisture - dry the skin surface by wiping away the moisture
- Color of the epidermis
- Skin texture
- Patient's skin condition is acceptable to continue therapy

Check the pad

Recheck the pad for the following at least every 30 minutes.

- position
- water status
- dry surface
- leaks
- cracks

Defibrillation considerations

1. Remove the pad to expose the patient's chest.
2. Remove excess moisture.

Disconnecting the pad from the hoses

1. Power off the controller.
2. Close pad clamps (Figure 4).
3. Unsnap and unlock the fitting. Pull the male coupling out to disconnect.
4. Discard pad after use per facility protocol.

Note: Do not discard the controller connector hose after use, save for next use.


Figure 4: Close the clamp

Figure 7: Snap

Figure 5: Insert

Figure 8: Lock

Figure 6: Twist

Figure 9: Open the clamp

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