Operations Manual Manual de uso Manuel d'utilisation

English (EN)
Español (ES)
Français (FR)

2022/01 1 8002-009-121 REV 1 www.C2Dx.com


Figure 1: Batch lot code
Figure 1: Código de lote
Figure 1: Code de lot


Figure 2: Connect one pad
Figura 2: Conecte una almohadilla
Figure 2 : Connecter un coussinet


Figure 3: Connect multiple pads
Figura 3: Conecte varias almohadillas
Figure 3 : Connecter plusieurs coussinets

Manufactured For:
C2Dx Inc.
555 E. Eliza St. | Ste. A
Schoolcraft, MI 49087
USA

C2Dx®

2022/01 8002-009-121 REV 1 www.C2Dx.com

EN 8002-009-121 REV 1

8002-009-121 REV 1 EN

Introduction

This manual assists you with the operation of the C2Dx Model 8002-062-222 T/Pad® hyper/hypothermia pad. Read this manual before operating this product and keep a copy on file. Set methods and procedures to educate and train your staff on the safe operation of this product.

Intended use

The T/Pump® localized therapy system supplies warm or cold water at controlled temperatures. Delivery of the water is through the thermal transfer devices (Multi-T-Pad®) for localized temperature therapy. The T/Pump is for use in situations where physicians determine heat therapy is necessary or desired.

Localized temperature therapy is of particular benefit in treating the following:

- Orthopedic conditions such as acute injuries, chronic pain, lower back, pain, muscle spasms, and strains;
- Skin trauma such as abrasions, boils, bruises, burns and contusions;
- Other medical conditions such as chronic arthritis, neuritis, phlebitis, tendonitis, and IV infiltration; and symptoms such as infections and localized pain.

Product description

The T/Pad is a heavy duty polymer. The Click-Tite® connectors provide the means to connect to the temperature controller. The button design allows water to flow when the pad is folded to form a customized fit.

Expected life

The T/Pad has a 90 day reusable expected life upon first use and under normal use conditions.

Contraindications

For heating:

- Application to a body surface with compromised blood flow (Ischemia area under pressure, arterial insufficiency).
- Application to a patient with an increased tendency to bleeding (aggravates potential for hemorrhage).
- Application to a body surface with possibility of malignancy (tissue metabolism is increased and therefore, the growth potential of the malignant tissue is increased).
- Treatment of hemotoma within first 24-48 hours (potential for re-bleeding and hemorrhage). Recent sprain or fracture (acute inflammatory response).
- In combination with topical solutions whose toxicity may be affected by the application of heat.
- In combination with other heat sources.

For cooling:

- Application to body surface with compromised blood flow (Ischemia, area under pressure, arterial insufficiency).
- Application to body surface with known vascular impairment such as frostbite, arteriosclerosis or ischemia
- Application to body surface in people with hypersensitivity to cold such as people with Raynaud's phenomenon, cold urticaria, cryoglobulinemia, and paroxysmal cold hemoglobinuria
- Application to body surface in people with impaired sensation
- In combination with topical solutions whose toxicity may be affected by the application of cold

Specifications

Material	rayon, etileno vinílico acetato (EVA), polipropileno, polietileno
Model 8002-062-222	
Length	22 in.
Width	15 in.
Compatible controllers T/Pump series temperature controllers	
Environmental Conditions	
Operating Temperature	60 °F (15,6 °C) - 90 °F (32,2 °C)
Relative humidity	30% - 75% - 95%
Storage and transportation	
Temperature	-20 °F (-29 °C) - 120 °F (48 °C)
Humidity relative	
30%	75%
10%	95%

C2Dx reserves the right to change specifications without notice.

Symbols

	General warning		Caution
	Consult instructions for use		Do not puncture
	Quantity		Single patient use
	Lot		Catalogue number
	Date of manufacture		Manufacturer

Contact information

Contact C2Dx Customer Service or Technical Support at 1-888-902-2239.

Warranty

C2Dx warrants that its models 8002-062-222 T/Pad products will be free from defects in manufacturing and workmanship for a maximum of 30 day use period. "Use" period begins when such devices individual package is opened. C2Dx's obligation under this warranty is expressly limited to supplying a product replacement, at its option, any product which is, at the sole discretion of C2Dx, found to be defective. If requested by C2Dx, products for which a claim is made shall be returned prepaid to the factory. Any improper use or alterations or repair by others in such manner as in C2Dx's judgement affects the product materially and adversely shall void this warranty. Any repair of C2Dx products using part not provided or authorized by C2Dx shall void this warranty. No employee or representative of C2Dx is authorized to change or waive any terms of this warranty.

Waiver and limitation

The express warranty set forth herein is the only warranty applicable to the product. Any and all other warranties, whether express or implied, including any implied warranty of merchantability or fitness for a particular purpose are expressly excluded by C2Dx. In no event shall C2Dx be liable for incidental or consequential damages.

Return policy

Products cannot be returned without prior approval from the C2Dx Customer Service Department. An authorization number will be provided which must be printed on the returned product. C2Dx reserves the right to charge shipping and restocking fees on returned product. Special, modified, or discontinued products are not subject to return.

Damaged merchandise

ICC regulations require that damage for damaged product must be made with written proof of receipt of the product. Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt. Upon prompt notification, C2Dx will file a freight claim with the appropriate carrier for damaged incurred. Claims will be limited to amount up to the actual replacement cost. In the event that this information is not available by C2Dx within the first 15 days of use following the delivery of the product, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full within thirty (30) days of receipt. Claims for any incomplete shipments must be made within thirty (30) days of invoice.

CAUTION

- U.S. Federal Law restricts this device to sale or by the order of a licensed physician.
- Do not use sharp objects or pins with this product.
- Do not clean this product if soiled; dispose of product according to hospital policy.
- Always inspect the product for tears, holes, stains, or any other damage before use. Always monitor more frequently for pediatric patients and with impaired circulation.
- Always use minimal layers of sheeting and incontinence pads. Too many layers between the patient's skin and the pad will reduce the cooling or warming capabilities of the system.

International warranty clause

This warranty affects U.S. domestic policy. Warranty outside the U.S. may vary by country. Please contact your local C2Dx Inc. representative for extra information.

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Assessing the patient's skin condition

Note: Follow your facility's skin care protocol. The following instructions are general guidelines.

- Document initial skin assessment before therapy.
- Check for any foreign items such as medical patches, IV's, or sutures.
- Check for fluids such as moisture or gels.
- Verify the fluid size is correct and pre-filled before applying to the patient.

Note: Apply the pad directly to the skin or add minimal layers.

Connecting and pre-filling the pad (connector hoses and pad)

1. Close the clamp on the connector hoses and pad (Figure 4).

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2. Check for any foreign items such as medical patches, IV's, or sutures.

3. Check for fluids such as moisture or gels.

4. Verify the fluid size is correct and pre-filled before applying to the patient.

5. Open all clamps on the connector hose and the pad (Figure 9).

6. Insert the power cord into a wall outlet.

7. Press the power button to turn on the controller. See the controller instructions for use.

8. Check that the pad fits well.

Notes:

- Iodine based disinfectants will stain the product. Staining does not affect the patient nor the use of the product.
- Do not use pins or sharp objects with this product.
- If you fold the pad, place the end of the pad that contains the tubing away from the patient.
- Interconnect using the Click-Tite® connectors to provide therapy to more than one body site at a time (Figure 2 and Figure 3).

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