
R_x ONLY

ENGLISH (EN)
ESPAÑOL (ES)
FRANÇAIS (FR)
PORTUGUÊS (PT)
日本語 (JA)

Contents

Introduction 2
Indications For Use 2
Contraindications 2
Description 2
For Use With 2
Safety Directives 3
Accessories 3
Definitions 3
Instructions 4
Quick Pressure Monitor Set 4
Indwelling Slit Catheter 5
Suggested Set-Up for the Pressure Monitor with Indwelling Slit Catheter Set 6
Troubleshooting 7
Cleaning 8
Recommended Equipment and Materials 8
To Wipe Down the Equipment 8
Testing and Maintenance 8
Battery Replacement 8
Calibration Verification 8
Storage and Handling 9
Disposal/Recycle 9
Specifications 10

Introduction

This instructions for use manual contains information intended to ensure the safe, effective, and compliant use of your product. This manual is intended for in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians. Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, including safety information, in-service training, or current literature, contact your C2Dx sales representative or call C2Dx customer service. Outside the US, contact your nearest C2Dx subsidiary.

Indications For Use

The C2Dx Intra-Compartmental Pressure Monitor System is intended to measure compartmental pressures to aid in the diagnosis of compartmental syndrome. Compartmental syndrome is defined as a condition in which the circulation and function of tissues within a closed space are compromised by increased pressure within the space. Excessive pressure must be relieved in order to prevent muscle and nerve necrosis, which lead to permanent loss of function and limb contracture.

Contraindications

None known.

Description

The Intra-Compartmental Pressure Monitor System is a hand-held, portable pressure monitor system specifically designed to measure tissue fluid pressures.

As an interpretive tool, the Intra-Compartmental Pressure Monitor System should be used along with other clinical indicators to aid the surgeon in the diagnosis of compartmental syndrome.

For Use With

This section identifies components intended to be used with the equipment to obtain a safe combination. The Intra-Compartmental Pressure Monitor System consists of three major components:

DESCRIPTION	REF
Pressure Monitor and Battery	0295-001-000
Quick Pressure Monitor Set <i>Sterile, disposable set includes a syringe pre-filled with saline, a side-ported needle, and a diaphragm chamber to maintain a sterile fluid pathway.</i>	0295-002-000
Indwelling Slit Catheter Set <i>Sterile, disposable set used for continuous pressure monitoring; set includes a slit catheter, a break-away needle, and extension tubing.</i> NOTE: The Quick Pressure Monitor Set (REF 0295-002-000) is required to use the Indwelling Slit Catheter Set.	0295-006-000

Safety Directives

WARNINGS:

- Before using this equipment, or any component compatible with this equipment, read and understand the instructions for use. Pay particular attention to safety information. Become familiar with the equipment before use.
- Only healthcare professionals trained and experienced in the use of this medical device should operate this equipment.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. C2Dx, as a manufacturer, does not recommend surgical procedure or technique.
- Upon initial receipt and before each use, inspect each component for damage. DO NOT use any equipment if damage is apparent.
- Upon initial receipt and before each use, clean the equipment as indicated. See the *Cleaning* section.
- The Pressure Monitor must remain dry. If liquid is spilled on the unit, the pressure reading may no longer be accurate. Discontinue use and return to C2Dx for service.
- DO NOT use this equipment in areas in which flammable anesthetics or flammable agents are mixed with air, oxygen or nitrous oxide.
- Take special precautions regarding electromagnetic compatibility (EMC) when using this medical electrical equipment. Install and place the equipment into service according to the EMC information contained in this manual. Portable and mobile radio frequency (RF) communications can affect the function of the equipment.
- ALWAYS operate the equipment within the specified environmental condition values. See the *Specifications* section. The reading of compartment pressure by the C2Dx Pressure Monitor System is relative to atmospheric pressure. It may vary versus an absolute pressure reading, which depends on the atmospheric pressure at the time the reading is taken. Readings should not be taken in environments with changing atmospheric pressure such as aircraft, hyperbaric chambers or any vehicle undergoing noticeable changes in elevation.
- Schedule to check the accuracy of the Pressure Monitor per hospital protocol or at least once a year. To test the accuracy, see the *Calibration Verification* section.

Accessories

WARNINGS:

- Use only C2Dx-approved accessories. Other accessories may result in increased emissions or decreased immunity of the system.
- DO NOT modify any equipment without the authorization of the manufacturer.
- DO NOT reuse, reprocess, or repack a device that is intended for single use only.
 - A single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing.
 - Design features may make cleaning difficult.
 - Reuse may create a contamination risk and compromise structural integrity resulting in operational failure.
 - Critical product information may be lost during repackaging. Failure to comply may lead to infection or cross infection and result in patient and/or healthcare staff injury.

This section identifies equipment that may be purchased to replace consumed, worn, or damaged equipment. Optional equipment may also be purchased to aid in operation.

The following C2Dx-approved accessories are sold separately:

DESCRIPTION	REF
Side-Ported Needle, 18 GA x 2.4 inch (6.1 cm)	0295-003-000
3 mL NaCl Fill/Syringe	0295-005-000
<i>Sterile, sodium chloride injection USP 0.9%. A vehicle, solvent, or diluent for parenterals.</i>	
Battery Door	0295-001-004
Housing Cover	0295-001-005
Storage Case	0295-008-000
9 V Alkaline Battery	0295-009-000

NOTE: For a complete list of accessories, contact your C2Dx sales representative or call C2Dx customer service. Outside the US, contact your nearest C2Dx subsidiary.

Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

SYMBOL	DEFINITION
	Power on
○	Power off
ZERO	Press this button to reset the indicator to a "00" value.
mmHg	Millimeters of mercury (unit of pressure measurement)

Instructions

Quick Pressure Monitor Set

WARNINGS:

- ALWAYS follow sterile procedures.
 - Before use of sterile accessories, inspect the package for damage and confirm the integrity of the sterile barrier. DO NOT use any accessory if damage is apparent or the sterile barrier has been compromised. The syringe solution and fluid path is sterile if the protective caps are in place and the package is intact.
 - Do not use the pre-filled syringe if the solution is discolored or contains a precipitate.
 - Keep the syringe at controlled room temperature. Keep from freezing.
 - The disposable syringe is for single use only. Discard any unused portion.
 - DO NOT resterilize or reuse the needles, syringe, or catheter. Discard these components after one use.
1. Turn the unit on. “_ _” is displayed for 5 seconds and then disappears. The device is now ready for use. If “_ _” is displayed continuously, see the *Troubleshooting* section.
 2. Remove the contents of the disposable pouch.
 3. Place the needle firmly on the tapered chamber stem.
 4. Remove the cap on the pre-filled syringe, and screw the syringe onto the remaining chamber stem. Be careful not to contaminate the fluid pathway. See Figure 1.


Figure 1

5. Open the monitor cover and inspect the transducer port for damage or contamination.
 - If the port is damaged, do not use the monitor. Return the device to C2Dx.
 - Remove any contamination. See the *Cleaning* section.
6. Place the chamber into the transducer well with the black surface facing down. Gently push until the disposable seats. See Figure 2.


Figure 2

7. Snap the cover closed — DO NOT FORCE. Make sure the latch has snapped.

WARNING: DO NOT use the monitor if the monitor cover is physically damaged or unable to snap closed.

NOTE: If the cover will not close, rotate the syringe so that a straight side of the flange rests against the base of the monitor. This will reduce the height of the flange to allow the cover to close.

8. Hold the needle approximately 45 degrees up from horizontal. Slowly force fluid through the disposable to purge it of air.

CAUTION: DO NOT allow saline to roll down the needle into the transducer well.

WARNING: To ensure accuracy, the device must be at the correct angle of insertion before the ZERO button is pressed.

9. Approximate the intended angle of insertion of the needle into the skin. To ensure accuracy, the display must read -8 to +9 mmHg. Press the ZERO button. The display will read “00” after approximately 2 seconds. See Figure 3.


Figure 3

NOTE: The display must read “00” before continuing. If the display does not read “00,” see the *Troubleshooting* section.

WARNING: DO NOT use the pressure readings when the display indicates low battery. Replace the battery (see the *Battery Replacement* section) and repeat the steps to obtain the pressure reading.

10. Insert the needle and slowly inject less than 3/10 mL of saline into the compartment to equilibrate with interstitial fluids.
11. Wait for the display to reach equilibrium.
12. Read the pressure.
13. For additional measurements, repeat steps 9 through 12. Make sure the unit is reset to zero.

Indwelling Slit Catheter

WARNING: The Indwelling Slit Catheter procedure must not exceed a 24-hour duration.

NOTE: ALWAYS start the procedure with a full syringe.

1. Perform steps 1 through 7 of the *Quick Pressure Monitor Set* instructions.
2. Remove the needle from the monitor.
3. Attach the slit catheter to the tapered chamber stem. Make sure the fit is very snug. See Figure 4.

NOTE: Optional extension tubing may be used between the slit catheter and the monitor if additional length is desired.


Figure 4

4. Remove the needle guard from the break-away needle, and push the catheter in until it is just visible at the end of the needle.

CAUTION: DO NOT let the catheter protrude from the needle, as it may be cut by the needle upon insertion.

5. Hold the needle approximately 45 degrees up from horizontal. Slowly force fluid through the disposable to purge it of air.
6. Bring the catheter tip (slit end) to the same height as the disposable diaphragm. See Figure 5. The display should read -8 to +9 mmHg.


Figure 5

WARNING: To ensure accuracy, the device must be at the correct angle of insertion (with the disposable diaphragm at the same height as the catheter tip) before the ZERO button is pressed.

7. Press the ZERO button.

NOTE: The display must read "00" before continuing. If the display does not read "00," see the *Troubleshooting* section.

WARNING: DO NOT use the pressure readings when the display indicates low battery. Replace the battery (see the *Battery Replacement* section) and repeat the steps to obtain the pressure reading.

8. Hold the catheter and the needle together and push into the compartment. Make sure the catheter's tip is protected inside the needle.
9. Push the catheter forward while pulling back on the needle. See Figure 6. Tape the catheter in place after the needle is out of the skin.


Figure 6

10. Remove the break-away needle.

NOTE: Pinch together the green fins on the needle, then pull the fins apart and split the needle away. See Figure 7.


Figure 7

11. Keep the chamber in the monitor at the same height as the estimated height of the slit tip of the catheter in the body. See *Suggested Set-Up for the Pressure Monitor with Indwelling Slit Catheter Set* section.

WARNING: The monitor must be securely taped to the patient to prevent its displacement and/or patient injury.


12. Inject less than 3/10 mL of saline into the compartment to equilibrate with interstitial fluids.
13. Wait for the display to reach equilibrium.
14. Record the pressure readings at the appropriate intervals. Repeat step 12 to ensure patency for subsequent measurements.

Suggested Set-Up for the Pressure Monitor with Indwelling Slit Catheter Set

Acute and Sports Medicine Set-Up

NOTE: Optional extension tubing not shown but may be necessary.


Patient shown standing.

Alternate Sports Medicine Set-Up Only

NOTE: Do not use this set-up for a patient lying down.


Patient shown standing.

Acute and Sports Medicine Set-Up

NOTE: Optional extension tubing not shown but may be necessary.


Patient shown lying down.

Alternate Acute Set-Up Only

NOTE: This set-up cannot be used for a patient standing up.

Patient shown lying down.

Troubleshooting

WARNING: DO NOT disassemble or service this equipment without the authorization of the manufacturer.

PROBLEM	CAUSE	ACTION
" _ _ " is displayed continuously.	The device requires repair.	Return the device to C2Dx.
"-E" is displayed.	The disposable diaphragm is higher than the needle or the catheter tip.	Correct the height condition.
	The device has been zeroed under pressure and the "-E" is still displayed when the pressure is removed.	Remove the disposable diaphragm (if present) and press the ZERO button. The device should read "00." Perform the procedure again with a new disposable diaphragm.
"-E" is displayed after pressing the ZERO button with no disposable diaphragm mounted.	The device requires repair.	Return the device to C2Dx.
Flashing "EO" is displayed.	The disposable diaphragm is not at the same height as the needle or the catheter tip during zeroing.	Correct the height condition then press the ZERO button.
	Fluid is present in the transducer port (the disposable diaphragm mounts to this metal fitting).	Allow the fluid to evaporate for 24 to 48 hours. CAUTION: Never insert any object into the transducer port.
	The disposable diaphragm is incorrectly purged before being mounted.	Perform the procedure again with a new disposable diaphragm.
	The ZERO button is accidentally pressed during measurement.	Turn the device off and on to restore the previous zero setting.
Flashing "EO" is displayed after pressing the ZERO button with no disposable diaphragm mounted.	The device requires repair.	Return the device to C2Dx.
More than 9 mmHg is displayed after the disposable diaphragm is mounted.	The needle or the catheter is blocked during purging.	Purge less than 3/10 mL of fluid through the disposable set and check for a blocked needle or catheter. Replace the disposable set if purging does not result in a reading of less than 9 mmHg.
	Fluid is present in the transducer port (the disposable diaphragm mounts to this metal fitting).	Allow the fluid to evaporate for 24 to 48 hours. CAUTION: Never insert any object into the transducer port.
Flashing "199" is displayed.	Pressure greater than 198 mmHg is present due to a blocked or kinked catheter or a blocked needle.	Check for a blocked or kinked catheter or a blocked needle. Replace the disposable set and perform the procedure again.
	The device requires repair.	Return the device to C2Dx.
"BAT" indicator is displayed.	The battery is depleted.	Replace the 9 V alkaline battery. See the <i>Battery Replacement</i> section.
Flashing "00," "60," or "EP" is displayed.	The device requires repair.	Return the device to C2Dx.
The monitor cover is loose or unable to snap closed.	The monitor cover is damaged.	Order a replacement Housing Cover (REF 0295-001-005). Do not use the monitor until the replacement cover is received and properly installed.
The device experiences sporadic electrical interference.	Electrical noise is present.	Turn off all electrical equipment not in use.
		Relocate electrical equipment and/or increase spatial distance between electrical equipment.
		Plug equipment into different outlets.

Cleaning

WARNINGS:

- ALWAYS clean the equipment as indicated upon initial receipt and before each use.
- DO NOT resterilize or reuse the needles, syringe, or catheter. Discard these components after one use.

CAUTIONS:

- DO NOT immerse any component in liquid. DO NOT allow liquids or moisture to enter any electrical connection.
- DO NOT sterilize the Pressure Monitor.
- DO NOT use solvents, lubricants, or other chemicals, including glutaraldehyde or similar chemical cleaners, unless otherwise specified.
- DO NOT use unapproved disinfectants. Failure to comply may cause equipment damage.

Recommended Equipment and Materials

- PPE as recommended by the disinfectant supplier
- Soft, lint-free cloth
- United States Environmental Protection Agency (US EPA) registered disinfectant with a claim for activity against Hepatitis B. The following disinfectants have been validated for use on the exterior surfaces of the equipment:

Sodium Hypochlorite Based - Clorox® Clean-Up® Disinfectant Cleaner with Bleach (EPA Reg. #67619-17)

Quaternary Ammonium Based - CaviCide® (EPA Reg. #46781-6)

To Wipe Down the Equipment

1. Turn off the monitor.
2. Wipe all external surfaces of the equipment with a soft, lint-free cloth moistened with a non-abrasive, hospital disinfectant prepared according to the manufacturer's instructions.
3. Upon the removal of visible gross soil, use a clean cloth moistened with disinfectant and wipe all surfaces. Surfaces must remain visibly wet at room temperature for at least the minimum time specified in the instructions for use supplied by the disinfectant manufacturer.
4. Remove any excess disinfectant solution using a soft, lint-free cloth moistened with water if required by the instructions supplied by the disinfectant manufacturer.
5. Inspect the equipment.

Testing and Maintenance

WARNINGS:

- Only individuals trained and experienced in the maintenance of this reusable medical device should test and maintain this equipment.
- DO NOT attempt to remove the screws on the bottom of the monitor and open the device. The electronics are coated for moisture resistance and are not serviceable.

NOTE: For service, contact your C2Dx sales representative or call C2Dx customer service. Outside the US, contact your nearest C2Dx subsidiary.

Battery Replacement

1. Turn off the monitor.
2. Remove the back cover of the monitor with a paper clip to expose the battery.
3. Pull the battery out of the unit.
4. Unsnap the battery from the connector.
5. Replace the battery with a new, 9 V alkaline battery.
6. Tuck the battery and the wires back inside the monitor.
7. Replace the cover, and snap the cover into place.

Calibration Verification

Verify the calibration of the Pressure Monitor (REF 0295-001-000) per hospital protocol or at least once per year.

CAUTION: Do not attempt to repair this unit. If recalibration is necessary, return the Pressure Monitor to C2Dx.

There are two methods of calibration verification: The Water Method (with a disposable unit), and the Air Method (without a disposable unit).

NOTE: The Water Method is preferred because it accounts for preloading characteristics of the monitor.

Water Method

Suggested Manometer Set-Up

NOTES:

- Use a solid water column manometer that measures in mmHg. An air-over-water manometer is not acceptable. Purge all the air bubbles from the system.
 - The Pressure Monitor should be checked with fluid in the disposable and not with air.
 - Fill the manometer with distilled water, which is free of contaminants. Do not use saline solution.
1. Install a fresh battery. See the *Battery Replacement* section.
 2. Prepare the Pressure Monitor with a disposable Quick Pressure Monitor Set. Follow steps 1 through 8 in the *Quick Pressure Monitor Set* section, but substitute distilled water for saline in the syringe.
 3. Set the Pressure Monitor flat on the table next to the manometer.
 4. Zero the unit, even if the display already reads "00."
 5. Insert the needle into the manometer so that the distal end of the needle is level with the zero point on the manometer.
 6. Compare the readings on the manometer and on the Pressure Monitor display at 25 mmHg, 50 mmHg, and 75 mmHg. The Pressure Monitor readings must be +/- 3.4 mmHg of the manometer readings.

NOTE: Keep the Pressure Monitor level. Tilting the unit can result in deviations in readings.

7. If the unit fails, repeat the procedure with another disposable Quick Pressure Monitor Set.
8. If recalibration is necessary, return the Pressure Monitor to C2Dx.

Air Method

1. Install a fresh battery. See the *Battery Replacement* section.
2. Turn the unit on. " _ _ " is displayed for five seconds and then disappears. The device is now ready for use. If " _ _ " is displayed continuously, see the *Troubleshooting* section.
3. Zero the unit, even if the display already reads "00."
4. Attach a 1/4 inch (6.4 mm) silicone tube to the transducer port under the clear lid of the monitor. Do not use a disposable diaphragm.

5. Connect the tubing to a calibrated electronic calibration unit or to a pressure calibrating unit.
6. The following readings should be taken:

CALIBRATION UNIT	PRESSURE MONITOR
25 mmHg	26 +/- 3 mmHg
50 mmHg	52 +/- 3 mmHg
75 mmHg	78 +/- 3 mmHg

7. If the unit fails, repeat the procedure.
8. If recalibration is necessary, return the Pressure Monitor to C2Dx.

Storage and Handling

CAUTION: ALWAYS store the equipment within the specified environmental condition values throughout its useful life. See the *Specifications* section.

Disposal/Recycle

WARNINGS:

- ALWAYS follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.
- ALWAYS follow the current local regulations governing the handling and disposal of sharps.
- ALWAYS follow the current local regulations governing biohazard waste to safely handle and dispose of surgical waste.

Specifications

Model:	REF 0295-001-000 Pressure Monitor and Battery
Dimensions:	117.5 mm x 65 mm x 36.74 mm
Mass:	150 grams (with battery)
Degree of Accuracy:	+/- 3.4 mmHg
Power Supply:	Internally Powered 9 V  (Direct Current) Operation: Standard 9 V  Alkaline Battery
Battery Life:	48 hours continuous operation minimum
Leakage Current:	Current: <= 100 µA Normal Conditions <= 300 µA Single Fault Conditions
Equipment Type:	 Type BF Applied Part
Mode of Operation:	Continuous Operation
Applied Part(s):	The entire Intra-Compartmental Pressure Monitor System device as defined by the manufacturer. See the <i>Product Safety Certification</i> standards listed below.
Ingress Protection:	IPX0 Ordinary Equipment
Product Safety Certification:	 CSA International C US International Electrotechnical Commission IEC 60601-1:1988, <i>Medical Electrical Equipment — Part 1: General Requirements for Safety - Second Edition</i> Amendment 1 (1991); Amendment 2 (1995); Corrigendum 1 (1995) Canadian Standards Association CAN/CSA-C22.2 No. 601.1-M90, <i>Medical Electrical Equipment — Part 1: General Requirements for Safety</i> CAN/CSA C22.2 No. 601.1S1-94-M94, Supplement No. 1-94 to 601.1-M90 Underwriters Laboratories UL 2601, <i>Medical Electrical Equipment</i>

Environmental Conditions:



Operation

Storage and Transportation



Temperature Limitation:

Humidity Limitation:

Atmospheric Pressure Limitation:

Guidance and manufacturer's declaration - electromagnetic emissions		
The Intra-Compartmental Pressure Monitor System (REF 0295-001-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the Intra-Compartmental Pressure Monitor System (REF 0295-001-000) should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Intra-Compartmental Pressure Monitor System (REF 0295-001-000) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Intra-Compartmental Pressure Monitor System (REF 0295-001-000) is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	n/a	
Voltage fluctuations/flicker emissions IEC 61000-3-3	n/a	

Guidance and manufacturer's declaration - electromagnetic immunity			
The Intra-Compartmental Pressure Monitor System (REF 0295-001-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the Intra-Compartmental Pressure Monitor System (REF 0295-001-000) should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	n/a n/a	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Intra-Compartmental Pressure Monitor System (REF 0295-001-000), including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.17 \sqrt{P}$ $d = 1.17 \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = 2.33 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m)</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  <p>(Non-ionizing electromagnetic radiation)</p> </div>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	

NOTE 1: At 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity			
The Intra-Compartmental Pressure Monitor System (REF 0295-001-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the Intra-Compartmental Pressure Monitor System (REF 0295-001-000) should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±2, 4, 6 kV contact ±2, 4, 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	n/a n/a	
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	n/a n/a	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	n/a n/a n/a n/a	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the alternating current mains voltage prior to application of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the Intra-Compartmental Pressure Monitor System (REF 0295-001-000)			
The Intra-Compartmental Pressure Monitor System (REF 0295-001-000) is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Intra-Compartmental Pressure Monitor System (REF 0295-001-000) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Intra-Compartmental Pressure Monitor System (REF 0295-001-000) as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17 \sqrt{P}$	80 MHz to 800 MHz $d = 1.17 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33 \sqrt{P}$
0.01	n/a	0.12	0.23
0.1	n/a	0.37	0.74
1	n/a	1.1	2.3
10	n/a	3.7	7.4
100	n/a	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Manufactured for C2Dx Inc