Biodesign®

OTOLOGIC REPAIR GRAFT

FP0108-01F

BIODESIGN® OTOLOGIC REPAIR GRAFT

PRODUCT DESCRIPTION

The Cook® Biodesign® Otolologic Repair Graft is a porous biomaterial composed of laminated extracellular collagen matrix derived from porcine small intestinal submucosa (SIS). SIS is obtained from the intestine using a process that retains the natural composition of matrix molecules such as collagen (Types I, III, VI), glycosaminoglycans (hyaluronic acid, chondroitin sulfate A and B, heparin, and heparin sulfate), proteoglycans, and fibronectin.^{1,2,3}

INTENDED USE

The Cook® Biodesign® Otolologic Repair Graft is intended for use as an implant material to aid in surgical repairs and as an adjunct to aid in the natural healing process in various otologic procedures, including but not limited to myringoplasty and tympanoplasty.

The device is supplied sterile and is intended for one-time use.

Rx ONLY This symbol means the following:

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

This product is intended for use by trained medical professionals.

OTOLOGIC REPAIR GRAFT This symbol means the following: Otolologic Repair Graft

CONTRAINDICATIONS

This device is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

GENERAL

• Users should be familiar with surgical technique for various otologic procedures, including but not limited to myringoplasty and tympanoplasty.

PRECAUTIONS

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions of the device.
- The device is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard device if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
- Avoid packing external canal with adherent dressings or applying excessive pressure in the ear canal.
- Please take care when opening tray packaging to ensure that device remains seated in the tray.

POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical device materials in otologic procedures:

- Abscess formation
- Allergic reaction
- Calcification
- Cholesteatoma
- Excessive redness, pain, swelling, or blistering
- Fever
- Infection
- Inflammation (initial application of surgical device materials may be associated with transient, mild, localized inflammation)
- Mastoiditis
- Migration
- Persistence of perforation
- Recurrence
- Retraction pockets
- Seroma
- Squamous cysts
- Thickening of the tympanic membrane

STORAGE

Store in a clean place at room temperature. Do not place in freezer. Avoid excessive heat. Keep dry.

STERILIZATION

This device has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

Required Materials

- A sterile dish (kidney dish or other bowl)
- Sterile forceps
- Sterile scissors
- Hydration fluid: sterile room temperature saline or sterile lactated Ringer's solution

NOTE: Always handle the device using aseptic technique. Minimize contact with latex gloves.

PREPARATION

1. Prepare the wound bed using standard surgical techniques. Ensure the area is free of exudate and devitalized tissue. An initial excision or debridement of the site may be necessary to ensure the wound edges contain viable tissue.
2. Control bleeding prior to applying the Biodesign® Otolologic Repair Graft.

SELECTION

3. Choose the appropriately sized device based on defect size, ensuring that the device covers the entire defect surface area and extends slightly beyond the wound margins.
4. Aseptically remove the device from the packaging.
5. Cut the device as needed with sterile scissors to cover the defect.

APPLICATION

6. The device may be placed in a dry state. Alternatively, if device hydration prior to application is preferred, hydrate the device in hydration fluid for no longer than 1 minute.
7. Verify that the device extends beyond the margins of the defect to obtain the appropriate overlap. Reposition the device as necessary.
8. Following placement, ensure the device is adequately hydrated.
9. Complete the standard surgical procedure and secure the device with preferred method of fixation (e.g. otologic packing or other appropriate method).
10. Discard any unused portions of the device according to institutional guidelines for medical waste.

POSTOPERATIVE CARE

- To provide the best environment for tissue integration, provide the patient with a list of standard post-procedure recommendations. The following should be considered:
 - Avoid exposing surgical site to direct moisture for 3 weeks unless approved by physician.
 - Avoid exposing surgical site to rapid changes in pressure for 2 to 3 weeks unless approved by physician.

ASSESSMENT

NOTE: If a gel forms on the wound surface, do not attempt to forcibly remove it. Successful absorption of Biodesign may form a caramel-colored or off-white gel. Do not remove this gel by debridement. This caramelization contains extracellular matrix (ECM), which continues to replace deficient and missing ECM in the wound.

11. As healing occurs, sections of Biodesign may gradually peel. Carefully remove any remaining loose product around the edge as needed.
12. Gently cleanse the wound surface with sterile saline; leave the ECM gel intact.
13. Carefully reassess the wound and record healing progression and other relevant information.

REFERENCES

1. Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med* 2007;18:537-543.
2. Hodde JP, Badylak SF, Brightman AO, Voytik-Harbin SL. Glycosaminoglycan content of small intestinal submucosa: A bioscaffold for tissue replacement. *Tissue Eng* 1996;2:209-217.
3. Data on File, Cook Biotech, Inc.

| | | | |
|--|------------------|---|--------------------------|
| LOT | Batch code | MR | Magnetic resonance safe |
| REF | Catalogue number | Manufacturer | |
| Consult instructions for use | | MD | Medical device |
| Contains biological material of animal origin | | Double sterile barrier system | |
| Country of manufacture | | Single sterile barrier system with protective packaging inside | |
| Do not resterilize | | Single sterile barrier system | |
| Do not re-use | | Sterilized using ethylene oxide | |
| Do not use if package is damaged and consult instructions for use | | UDI | Unique device identifier |
| Keep dry | | Use-by date | |

MANUFACTURER
COOK BIOTECH
INCORPORATED
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.

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© COOK BIOTECH 2023
- Données internes chez Cook Biotech, Inc.