COOK**MEDICAL****CE**
0123

EN	Fixed Core Wire Guides
3	Instructions for Use
CS	Vodící dráty s pevným jádrem
4	Návod k použití
DA	Kateterledere med fast kerne
6	Brugsanvisning
DE	Führungsdrähte mit fester Seele
7	Gebrauchsanweisung
EL	Συρμάτινοι οδηγοί με σταθερό πυρήνα
9	Οδηγίες χρήσης
ES	Guías de alma fija
11	Instrucciones de uso
FR	Guides à âme centrale
12	Mode d'emploi
HU	Rögzített magú vezetődrotok
14	Használati utasítás
IT	Guide ad anima fissa
15	Istruzioni per l'uso
NL	Voerdraden met een vaste kern
17	Gebruiksaanwijzing
NO	Ledevaiere med fast kjerne
19	Bruksanvisning
PL	Prowadniki z nieruchomym rdzeniem
20	Instrukcja użycia
PT	Fios guia com parte central fixa
22	Instruções de utilização
SV	Ledare med fast kärna
23	Bruksanvisning

T _ F C W G _ R E V Z

FIXED CORE WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Fixed Core Wire Guides are made with a stainless steel inner mandril and stainless steel outer coil. They are available in a variety of diameters, lengths, and tip configurations. Some devices are provided with a green PTFE coating to ease coaxial use with other devices. Refer to the product label for product specifications.

INTENDED USE

Fixed Core Wire Guides are intended to facilitate the placement of devices during diagnostic and interventional procedures.

CONTRAINDICATIONS

None known

WARNINGS

- This product is not intended for use in the coronary arteries.
- This product is a delicate instrument. Avoid forceful angulation.
- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape or shear material from the wire guide.
- Altering the tip's configuration or curve manually may damage the wire guide.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of devices in angiographic procedures should be employed.
- Do not attempt to torque the wire guide.
- Use medical imaging when you manipulate the wire guide. Do not advance or manipulate the wire guide without visual evidence of the corresponding movement of the distal tip.
- When you use the wire guide with another device, consider the end-hole size and the length of the device in order to ensure a proper fit between the wire guide and the device.
- Do not advance or withdraw a wire guide when resistance is encountered, as a perforation could occur.
- Use only the flexible end as the initial inserted end. Inserting from the rigid end may cause damage to tissue or the device.
- Inspect the wire guide for kinks or damage prior to use.

INSTRUCTIONS FOR USE

1. Flush the wire guide holder by attaching a syringe with heparinized saline or sterile water to the fitting of the wire guide holder. Inject enough solution to wet the wire guide surface entirely.

NOTE: If flushing through the wire guide holder is not possible, remove the wire guide from the holder and place it in a bowl of heparinized saline or sterile water, or wet the wire guide surface over the entire length using gauze moistened with heparinized saline solution.

2. Carefully remove the wire guide from the holder.
3. If needed, insert a wire guide insertion tool (provided) through the valve assembly or hub of the guiding sheath or other interventional device. Insert the tip of the wire guide through the insertion tool.
4. Standard wire guide techniques may now be employed.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for more information.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd