EN  
3

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Bruksanvisning

T \_ S A L \_ R E V 3

## ADVANCE® SALIVARY BALLOON CATHETER WITH WIRE GUIDE

**CAUTION: U.S. federal law restricts this device to sale on or by the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

The Advance Salivary Balloon Catheter with Wire Guide consists of a balloon catheter and a wire guide.

The balloon catheter is a double-lumen catheter with a balloon near its distal tip. The catheter consists of two coaxial lumens, which are labeled "DISTAL" and "BALLOON." The distal lumen extends the length of the catheter and is used to advance the catheter over the supplied wire guide. The balloon lumen is used to expand the balloon. The balloon diameter (mm) and balloon length (mm) are inscribed on the tip of the manifold. Refer to product label for product specifications.

The balloon is manufactured from an extra-thin wall, high-strength, minimally compliant material. **Particular care should be taken in handling the balloon to prevent damage. It will inflate to the indicated size parameters when utilizing proper pressure recommendations. Adhere to balloon pressure parameters indicated. (See Compliance Card Insert.)** Use of a pressure gauge is recommended to monitor inflation pressures.

The wire guide is a stainless steel mandril with a stainless steel coil at the distal tip. The wire guide is 0.014 inch (0.36 mm) in diameter and 50 cm in length.

### INTENDED USE

The Advance Salivary Balloon Catheter is intended for the dilation of lesions and strictures in the submandibular and parotid salivary ducts.

### CONTRAINDICATIONS

- Acute glandular infection
- Acute purulent sialadenitis

### WARNINGS

- This product is not intended for use in the vasculature.
- This product is not intended for use in the sublingual salivary duct or glands.
- This product is not intended for use through the working channel of a sialendoscope, as this could result in shearing of the balloon.
- Do not exceed rated burst pressure.<sup>1</sup> Rupture of balloon may occur. **Adhere to balloon inflation pressure parameters. (See Compliance Card Insert.)** Overinflation may cause rupture of the balloon and may result in damage to the duct wall. Use of a pressure gauge is recommended to monitor inflation pressures.
- Do not preinflate the balloon before inserting it into the salivary duct. Doing so may make balloon introduction more difficult.
- Do not use a power injector for balloon inflation. Balloon rupture may occur.
- Do not use a power injector for injection of contrast medium through catheter lumen marked "DISTAL." Catheter rupture may occur.

<sup>1</sup>The burst pressure data were analyzed using factors for a one-sided tolerance to determine with 95% confidence that 99.9% of these balloons would not burst at or below calculated rated burst pressure.

### PRECAUTIONS

- The product is intended for use by physicians trained and experienced in sialendoscopy.
- When inserting, manipulating or withdrawing the device, do not use excessive force.
- Flush the distal catheter lumen if resistance is felt when passing the catheter over the wire guide.
- The balloon is constructed of a heat-sensitive material. Do not heat or attempt to shape the catheter tip.
- Use only saline or sterile water as inflation medium. Never use air or any gaseous medium to inflate the balloon.

- Use the device prior to the “Do Not Use After” date specified on the package.
- The device is not intended for the delivery of stents.
- During device withdrawal, take care to keep the shaft as straight as possible. Withdrawing the device at an excessive angle or curvature may damage the shaft.
- Do not attempt dilation if the wire guide and balloon cannot cross the lesion or stricture.
- Do not attempt usage of this product in the presence or suspected presence of unfavorable patient anatomy.

## POTENTIAL ADVERSE EVENTS

The following complications are possible with the use of balloon catheter devices:

- Pain
- Infection
- Bleeding
- Tissue dissection, perforation, rupture or injury
- Duct occlusion
- Tissue avulsion
- Damage to the facial or lingual nerve

## INSTRUCTIONS FOR USE

The Advance Salivary Balloon Catheter requires imaging in order to ensure proper placement and inflation of the balloon catheter. The choice of imaging depends on physician preference. **Table 1** shows the recommended imaging modalities and compatible ancillary devices.

**Table 1: Recommended Imaging Modality and Compatible Ancillary Devices**

Imaging Modality	Treatment Site	Introducer Sheath
Sialendoscopic Visualization*	Submandibular	6 Fr
Sialendoscopic Visualization with Ultrasound Imaging	Submandibular	6 Fr*
	Parotid	N/A

\*1.3 mm OD sialendoscope maximum

\*Alternatively, the balloon catheter and wire can be used without an introducer sheath.

## Balloon Catheter Preparation

1. Choose a balloon catheter diameter appropriate for treatment of the salivary duct. (See **Table 2.**)

**Table 2: Recommended Duct Sizing**

Reference Balloon Diameter	Recommended Minimum Duct Diameter	Recommended Minimum Duct Length
1.5 mm	1.5 mm	2 cm
2.0 mm	2.0 mm	
2.5 mm	2.5 mm	
3.0 mm	3.0 mm	

2. Remove the balloon catheter from the package. Inspect the balloon catheter to ensure that no damage has occurred during shipping.
3. Remove the yellow-tipped shipping stylet from the distal tip of the balloon catheter.
4. Remove the protective orange balloon sleeve from the balloon catheter.
5. Apply negative pressure to the lumen labeled "BALLOON" before introducing it into the salivary duct.

### **Balloon Catheter Introduction and Inflation with Introducer Sheath and Sialendoscopic Visualization**

#### **(For Use in the Submandibular Salivary Duct)**

1. Use an appropriately sized introducer sheath to access the desired salivary duct. (See **Table 1.**)
2. Use the sialendoscope to view the desired treatment site, and then pre-position the wire guide, floppy end first. If using a sialendoscope with an instrument channel, the wire guide may be inserted through the instrument channel. If there is no instrument channel, insert the wire guide through the introducer sheath into the salivary duct alongside the sialendoscope.
3. Advance the wire guide until it has passed through the desired treatment site.
4. While maintaining position of the wire guide within the salivary duct, remove the sialendoscope from the salivary duct and from the wire guide.
5. Place the balloon catheter over the wire guide. Advance the balloon catheter over the wire guide until the proximal end of the balloon catheter enters the introducer sheath.  
**NOTE:** In order to prevent trauma to the duct, be sure to maintain the position of the wire guide when you advance and manipulate the balloon catheter.
6. Insert the sialendoscope into the introducer sheath behind the balloon catheter and advance the sialendoscope through the introducer lumen.
7. Advance the balloon catheter and sialendoscope together until the treatment site is visualized.
8. Advance the balloon catheter independently until the balloon is positioned across the desired treatment site. **Do not advance the balloon into the smaller branch of duct as undesired dilation of the smaller branch may occur.**  
**NOTE:** If you meet with resistance while advancing the balloon catheter, do not force the balloon catheter.
9. To ensure that the sialendoscope is not alongside the balloon during inflation, maintain position of the balloon catheter while withdrawing the sialendoscope until the balloon catheter shaft is visible.
10. Prepare an inflation device with saline or sterile water and attach it to the balloon lumen.
11. Inflate the balloon to the desired pressure. **Adhere to recommended balloon inflation pressures. (See Compliance Card Insert)**
12. If the balloon loses pressure and/or ruptures, pull negative pressure on the balloon inflation port and remove the balloon and sialendoscope as a unit.

#### **Balloon Deflation and Withdrawal**

1. Use an inflation device or syringe to completely deflate the balloon. Allow adequate time for the balloon to deflate.
2. Maintain vacuum on the balloon and withdraw the balloon catheter.  
**NOTE:** Do not use excessive force when removing the balloon.  
**NOTE:** If you meet with resistance while using a sialendoscope, remove the sialendoscope first, and then the balloon catheter.
3. Confirm adequate treatment with sialendoscope after balloon removal.
4. Remove sialendoscope and wire guide.

## **Balloon Catheter Introduction and Inflation with Sialendoscopic Visualization and Ultrasound Imaging**

### **(For Use in the Submandibular or Parotid Salivary Ducts)**

**NOTE:** If targeting the submandibular salivary duct, an appropriately sized introducer sheath may be used if desired. (See **Table 1**)

1. Using a sialendoscope to visualize the desired treatment site, pre-position the wire guide, floppy end first. If using a sialendoscope with an instrument channel, the wire guide may be inserted through the instrument channel. If there is no instrument channel, insert the wire guide through the introducer sheath into the salivary duct alongside the sialendoscope.

**NOTE:** If an introducer sheath is not utilized, the papilla may be dilated before introducing the balloon catheter.

2. Advance the wire guide until it has passed through the desired treatment site.
3. Remove the sialendoscope, leaving the wire guide positioned through the desired treatment site.
4. Place the balloon catheter over the wire guide. While using ultrasound guidance, advance the balloon catheter over the wire guide until the balloon is positioned across the desired treatment site.

**NOTE:** In order to prevent trauma to the duct, be sure to maintain the position of the wire guide when you advance and manipulate the balloon catheter.

**NOTE:** If you meet with resistance while advancing the balloon catheter, do not force the balloon catheter.

**NOTE:** If necessary to better visualize balloon placement with ultrasound guidance, you may preinflate the balloon using low pressure (< 2 ATM) and deflate the balloon while it is in the duct. **Ensure the tip of the balloon is not positioned into the smaller branch of duct as undesired dilation of the smaller branch may occur.**

5. Prepare an inflation device with saline or sterile water and attach it to the balloon lumen.
6. Inflate the balloon to the desired pressure. **Adhere to recommended balloon inflation pressures. (See Compliance Card Insert)**
7. If the balloon loses pressure and/or ruptures, pull negative pressure on the balloon inflation port and remove the balloon catheter.

### **Balloon Deflation and Withdrawal**

1. Use an inflation device or syringe to completely deflate the balloon. Allow adequate time for the balloon to deflate.
2. Maintain vacuum on the balloon and withdraw the balloon catheter.  
**NOTE:** Do not use excessive force when removing balloon.
3. Confirm adequate treatment with sialendoscope after balloon removal.
4. Remove sialendoscope and wire guide.

### **HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

### **REFERENCES**

These instructions for use are based on experience from physicians. Refer to your local Cook sales representative for more information.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd