**COOK**  
MEDICAL

**CE**  
0123

**EN**  
3

**SialoCath™ Salivary Duct Catheter**

Instructions for Use

**CS**  
4

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T \_ S D C \_ R E V 2

# SIALOCATH™ SALIVARY DUCT CATHETER

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

## DEVICE DESCRIPTION

The SialoCath™ Salivary Duct Catheter is a single lumen catheter designed to assist in the treatment of salivary duct conditions.

## INTENDED USE

The SialoCath™ Salivary Duct Catheter is intended to aid in the removal and treatment of salivary duct stones and other duct conditions through flushing and/or fluid irrigation of the salivary ducts.

## CONTRAINDICATIONS

None known

## WARNINGS

- The product is not intended for use in the vasculature.
- This product is not intended for use in the sublingual salivary ducts or glands.
- This product is not intended for use through the working channel of a sialendoscope, as this could result in damage to the device.
- This product is not intended to be used with a power injection device.

## PRECAUTIONS

- The product is intended for use by physicians trained and experienced in sialendoscopy. Standard sialendoscopic techniques should be employed.
- Ensure compatibility of instruments and devices to be used with the SialoCath™ Salivary Duct Catheter.
- Do not attempt to insert or withdraw the catheter if resistance is felt.
- Do not attempt to insert the catheter in the presence or suspected presence of unfavorable patient anatomy.

## POTENTIAL ADVERSE EVENTS

Potential adverse events associated with the insertion and use of the SialoCath™ Salivary Duct Catheter include but are not limited to:

- Pain
- Infection
- Bleeding
- Tissue dissection, perforation, rupture or injury
- Duct occlusion
- Tissue avulsion
- Damage to facial or lingual nerve

## INSTRUCTIONS FOR USE

1. Upon removal of the device from the package, visually inspect with particular attention to kinks, bends, or breaks. If an abnormality that would restrict proper working conditions is apparent, do not use. Ensure the device is compatible with any instruments and/or devices desired to be used during the procedure.

2. Insert the catheter into the salivary duct to be treated over a 0.018" (0.46 mm) maximum diameter wire guide or through an appropriately sized introducer sheath.
3. Attach a syringe for flushing and/or fluid irrigation.  
**NOTE:** The catheter may be used after stone fragmentation to aid in the removal of stone segments.
4. After completion of flushing and/or fluid irrigation, remove the catheter from the salivary duct.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for more information.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd