

THE NEXT GENERATION HANDLE AND BLADES

Optimized for Your
Surgical Performance

NEW DISPOSABLE HANDLE

- New magnetic blade insertion
- New blade hub eliminates contact with blade
- Distinct Cut and Coag buttons
- Single-use ensures reliability

PRODUCT ORDERING INFORMATION

70139050	Scalpel Handle (quantity 1)
70239050	Scalpel Handle (quantity 6)
70235810	#10 Scalpel Blade (10 per box)
70135810	#10 Scalpel Blade (24 per box)
70235812	#12 Scalpel Blade (10 per box)
70235815	#15 Scalpel Blade (10 per box)
70135815	#15 Scalpel Blade (24 per box)

NEW BLADES

- Electro-chemically honed and 50% sharper than previous generation blades
- New non-stick coating virtually eliminates coagulum build up
- New blade hub insulates blade shaft, reducing risk of thermal injury

