



STIC INTRA-COMPARTMENTAL PRESSURE MONITOR

DISPOSABLES OFFERING



20-Gauge Side Ported Needle



18-Gauge Needle



Indwelling Slit Catheter Needle
for Continuous Monitoring

PRODUCT ORDERING INFORMATION

295-002-000

Quick STIC 18-gauge Needle Pressure Monitor Set (6/pkg).

Each set includes pre-filled syringe, diaphragm chamber and side-ported 18-gauge needle.

295-004-000

Quick STIC 20-gauge Needle Pressure Monitor Set (6/pkg).

Each set includes pre-filled syringe, diaphragm chamber and side-ported 20-gauge needle.

295-007-000

Side-ported, 20-gauge Needle (12/pkg).

295-006-000

Indwell Slit Catheter Set (12/pkg).

Each set includes slit catheter, breakaway needle and extension tubing.