

COMPARTMENT SYNDROME

A Timely Diagnosis is Critical




STIC PRESSURE MONITOR

OBJECTIVE INPUT FOR PROMPT, DATA-DRIVEN, MORE INFORMED DECISIONS

The **STIC Intra-Compartmental Pressure Monitor** provides quick and continuous readings of compartment pressure, adding valuable data to your clinical assessment. Formerly made by Stryker and now through C2Dx, STIC has been long recognized as the **gold standard** with more than **30 years** of clinical history.

THE OUTCOMES FOCUSED ON TIME



“As one of the major orthopedic emergencies, acute compartment syndrome carries significant risks when the diagnosis is delayed or missed.”

DELAYED TREATMENT

3x
MORTALITY &
AMPUTATION

↑ INFECTION
46%

4+
hrs ISCHEMIC
INJURY TO
MUSCLE/NERVES

EARLY DETECTION

5.5 ↓ DAYS
Overall

2.4 ↓ DAYS
in the ICU




* Data based on several clinical studies assessing how time to diagnosis and treatment affects patient outcomes.

Scan for video on how to use the device

The measurement of intra-compartmental pressure using the Gold Standard is the most reliable test for proper diagnosis.”2*

— Brit Long, MD

OVER 30 YEARS OF CLINICAL HISTORY



Quick and Simple Hand-Held Monitor for Accurate Measurement of Compartment Pressure



FAST: Simple, sterile, disposable



ACCURATE: Side-ported needle for proven accuracy and over 30 years of clinical history



VERSATILE: Provides both quick and continuous (up to 24hrs) pressure monitoring options



CONVENIENT: Pre-filled syringe, hand-held and easily transported



COST-EFFECTIVE: Low initial purchase and economical use. Compartment Syndrome diagnosis is a **reimbursable procedure - CPT code 20950**.

DOC-514-Rev 2

Copyright© 2021 C2Dx
Printed in U.S.A.

555 East Eliza St, Ste. A
Schoolcraft, MI 49087 USA
t: 888 902 2239

www.c2dx.com

PRODUCT ORDERING INFORMATION

STIC Pressure Monitor System

(includes the following-no substitutions or deletions allowed)

- 295-001-000 STIC Pressure Monitor
- 295-002-000 Quick Pressure Monitor Set (6/pkg). Each set includes pre-filled syringe, diaphragm chamber, and side-ported needle
- 295-006-000 Indwell Slit Catheter Set (12/pkg). Each set includes slit catheter, breakaway needle, and extension tubing.
- 295-009-000 9V Alkaline Battery
- 295-008-000 Storage Case

295-000-000

295-001-000 STIC Pressure Monitor

Quick STIC 18-gauge Needle Pressure Monitor Set (6/pkg)

295-002-000

Each set includes pre-filled syringe, diaphragm chamber and side-ported 18-gauge needle

Quick STIC 20-gauge Needle Pressure Monitor Set (6/pkg)

295-004-000

Each set includes pre-filled syringe, diaphragm chamber and side-ported 20-gauge needle

Indwell Slit Catheter Set (12/pkg)

295-006-000

Each set includes slit catheter, breakaway needle and extension tubing

ADDITIONAL ITEMS

295-007-000 Side-ported, 20-gauge Needle (12/pkg)

295-001-004 Battery Door

295-001-005 Housing Cover

295-009-000 9V Alkaline Battery

295-008-000 Storage Case

295-016-000 Pressure Monitor Wall Mount storage Cabinet

1. Little M.T.M., Lin C.A., Vrahas M.S. (2019) Legal Aspects of Compartment Syndrome.

In: Mauffrey C., Hak D., Martin III M. (eds) Compartment Syndrome. Springer, Cham

2. Long, Brit, MD, and Michael Gottlieb, MD, RDMS. The Dreaded Acute Compartment Syndrome. EmDOCs.net - Emergency Medicine Education. N.p., 02 Jan. 2019. Web.

*Data sources and the full list of references available upon request