

Extract salivary stones.




Illustration by Lisa Clark

NCircle[®]

TIPLESS SALIVARY STONE EXTRACTOR

Atraumatic design


The tipless design minimizes trauma to the duct and reduces the risk of ductal wall perforation.

Ergonomic

The handle is easy to use and allows you to deploy the stone extractor single-handedly.

Durable

The stone extractor's nitinol material is flexible and kink resistant.




NCircle®

TIPLESS SALIVARY STONE EXTRACTOR

Used for manipulation, entrapment and extraction of calculi and other foreign bodies in the salivary ducts, and to minimize migration of salivary stones during interventional procedures.

| Order Number | Reference Part Number | Catheter Fr | Catheter Length cm | Basket Diameter mm |
|--------------|-----------------------|-------------|--------------------|--------------------|
| G23797 | SSEB-1.5-115-10 | 1.5 | 115 | 10 |


Related products



Kolenda Salivary Access Introducer Set

Used for the introduction of diagnostic or interventional devices in the submandibular salivary duct.

| Order Number | Reference Part Number | Introducer ID mm/inch | Wire Guide Diameter inch | Sheath Length cm |
|--------------|-----------------------|-----------------------|--------------------------|------------------|
| G24129 | SDAS-6.0-18 | 2.2/.087 | .018 | 5.5 |
| G24132 | SDAS-5.0-18 | 1.9/.074 | .018 | 5.5 |



Wire guides

Used to gain and maintain access to the salivary duct.

| Order Number | Reference Part Number | Wire Guide Diameter inch | Wire Guide Length cm | Tip Configuration |
|--------------|-----------------------|--------------------------|----------------------|-------------------|
| G02274 | C-SF-15-50 | .015 | 50 | Straight |
| G02255 | C-SF-18-40 | .018 | 40 | Straight |

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4