


Hercules® 100

Transnasal Esophageal Balloon

Created for patients. Designed for ENTs.



MANUFACTURER

Cook Medical

Since 1963, Cook Medical has worked with physicians to develop minimally invasive technologies. Today, Cook offers medical devices, biologic materials, and cellular therapies to deliver better patient outcomes more efficiently. We remain family owned so we can focus on what we care about most: patients, our employees, and our communities.

PRODUCT DESCRIPTION

The Hercules 100 Transnasal Esophageal Balloon is a three-stage dilation balloon with a 100 cm catheter length that is comparable to most transnasal endoscopes. It is the first balloon designed for transnasal esophageal dilations. These dilations can be performed as in-office treatments.

INTENDED USE

The Hercules 100 Transnasal Esophageal Balloon is used to dilate strictures in the esophagus and the upper esophageal sphincter. The device is supplied sterile and is intended for single use. The Cook Sphere® Inflation Device is designed to be used with Hercules balloons.

DOCUMENTS

- [FDA clearance letter](#)
- [Instructions for Use](#)
- [Latex letter](#)
- [Reimbursement codes](#)

MARKETING MATERIALS

- [Product brochure](#)
- [Product datasheet](#)
- [Bibliography](#)
- [Procedural animation](#)
- [Procedural video](#)

PRODUCT WEBLINK

[CookMedical.com/Laryngology](https://www.cookmedical.com/Laryngology)

For more information,
please contact your local
Cook sales representative.

