

Dilate salivary duct strictures.




Illustration by Lisa Clark

Advance[®]

SALIVARY DUCT BALLOON CATHETER SET

High-strength

The high-strength, minimally compliant material allows for the complete dilation of the stricture.

Minimally invasive

The flexible, low-profile design supports minimally invasive intervention through radial dilation.

Multi-modal visualization

The catheter is compatible with direct endoscopic visualization or ultrasound imaging based on physician preference.

C2Dx[®]

c2dx.com



Advance[®]

SALIVARY DUCT BALLOON CATHETER SET

Used for the dilation of strictures in the submandibular and parotid salivary ducts.

The balloon is constructed of an extra-thin wall and high-strength material that is intended to be used with sialendoscopic visualization in the submandibular salivary ducts.

Components

- Advance Salivary Duct Balloon Catheter
- Wire guide

Order Number	Reference Part Number	Wire Guide Diameter inch	Inflated Balloon Diameter mm	Balloon Length cm	Preinflated Balloon Diameter mm
G34991	SDBS-14-30-1.5-2	.014	1.5	2	1
G34992	SDBS-14-30-2.0-2	.014	2.0	2	1
G34993	SDBS-14-30-2.5-2	.014	2.5	2	1
G34994	SDBS-14-30-3.0-2	.014	3.0	2	1

Related products



Sphere[®] Inflation Device

Order Number	Reference Part Number	Fluid Capacity mL	Maximum Pressure atm
G34903	CID-60-15	60	15

Some products or part numbers may not be available in all markets. Contact your local C2Dx sales representative or Customer Service for details.



* Copyright© 2024 C2Dx
 * Printed in U.S.A.
 * 555 East Eliza St, Ste. A
 Schoolcraft, MI 49087 USA
 * t: 888 902 2239
 www.c2dx.com